

FORM NRH 21B Effective Date: November 28, 2016

NEBRASKA DEPARTMENT OF HEALTH AND HUMAN SERVICES
DIVISION OF PUBLIC HEALTH
X-RAY PROGRAM

DENTAL EQUIPMENT PERFORMANCE EVALUATION

Evaluation Date:					
Registrant:		Re	egistration Number:		
Service Company:	Re	Registration Number: RS			
Survey Instrument:		Ca	Calibration Date:		
Ion Chamber: 🗌 Wit	hin A Housing	External Probe			
Control Panel Informat	ion				
Manufacturer:		Model Nur	nber:		
Serial Number:	<u>.</u>	Location:			
TIMER The accuracy of the time not obtainable, the time	er accuracy must be ±	10% of the indicated tir	ne with the testing per		
Manufacturer's Specific					
Time Used for Testing:	·		s 🗌 milliseconds 🗌 pı	ulses	
Measured:					
Deviation:	%	%	%	%	
EXPOSURE REPR When all technique factors automatic exposure contechniques. $C = \frac{s}{\overline{x}} = \frac{1}{\overline{x}} \left[\frac{\sum_{i=1}^{n} (x_i)}{n-1} \right]$	$\left \frac{-\overline{x}}{2} \right ^{2} = \left \frac{1}{2} \right ^{1/2}$ s = es X = m X = m X = ith		uirement applies to cli tion of the population ns in sample		
Technical Factors Sele	ected: kVp _	mA / mAs	seconds	/ milliseconds / pulses	
Measured:	mR	mR	mR	mR	
Coefficient of Var Manufacturer:		Serial Nur	nber:		



FORM NRH 21B Effective Date: November 28, 2016

KILOVOLT PEAK

PASS

FAIL

If the registrant possesses documentation of the appropriate manufacturer's kilovolt peak specifications, the radiation machine must meet those specifications. If the registrant does not possess documentation of the appropriate manufacturer's kilovolt peak specifications, the indicated kilovolt peak must be accurate to within ±10% of the indicated setting(s). For dental radiation generating equipment with fewer than three fixed kilovolt peak settings, the radiation machine will be checked at those settings.

Manufacturer's Specifications: _____ OR ±10%

Kilovolt Peak Used for Testing:

	-			
Measured:				
Deviation:	%	%	%	%
TUBE STABILITY PASS FAIL The tube must remain physically stable during exposures. In cases where tubes are designed to move during exposure, the registrant must assure proper and free movement of the unit. FAIL				
COLLIMATION Field limitation must m	eet the requirements c	of 180 NAC 21-007.07	PASS or 21-007.08	
Intraoral: Source-to-Skin Distand <i>AND</i> Beam limited to a diam	. ,		an 18 cm:	🗌 No
Panoramic: Transverse axis: x-ray <i>AND</i> Vertical axis: x-ray bea				🗌 Yes 🗌 No
Cephalometric: Rectangular collimation: x-ray field does not exceed 2.0% of the source-to-image receptor distance for the length or the width of the image receptor: Yes No OR Circular or polygon collimation: x-ray field does not exceed 2.0% of the source-to-image distance for the				
diagonal of the image receptor: <pre> Yes</pre> No IN-AIR EXPOSURE MEASUREMENT				
For use on dental intra	oral systems on an av	erage adult patient thic	kness in routine bitewi	ng radiography
Technique Factors:	kVp	mA / mAs	seconds / milli	iseconds / pulses
For Calculated Measure Source to Skin Distance		cm Source to	Detector Distance (SD	D):cm
Measured:	mR 🔤 0	Calculated Measureme	nt 🗌 Direc	t Measurement

Surveyor Name:	Surveyor Signature:	
carveyer name.	 ourreyer eignature.	



DETERMINING IN-AIR EXPOSURE MEASUREMENT FOR INTRAORAL DENTAL EXAMINATIONS

A. CALCULATION

Note: Ion chambers may be located within the instrument housing rather than within an external probe. In this situation, the distance from the top surface of the housing to the ion chamber below must be known. If this type of instrument is used for the measurements, the inverse square law must be utilized for accurate results.

IAE = Measured X (SDD ÷ SSD)²

Where:

IAE = in-air exposure Measured = indicated exposure on measuring instrument SDD = source (target) to detector (ion chamber) distance in centimeters SSD = source (target) to skin distance in centimeters

- 1. Place the tip of the cone within $\frac{1}{2}$ inch from the housing of the measuring instrument.
- 2. Measure the distance from the source to the entrance/tube side surface of the housing. This is the SSD.
- 3. Determine the distance from the source to the ion chamber within the housing. This is the SDD.
- 4. Select the kVp, mA(s), and time normally used for an average adult patient thickness in routine bitewing radiography at that facility.
- 5. Make an exposure and document the radiation output in millirem.
- 6. Using the above formula, calculate the in-air exposure.

B. DIRECT MEASUREMENT

Note: Use this procedure only if an external probe (ion chamber) is available for the measurements.

- 1. Position the tube so the end of the cone is not greater than ½ inch from the probe. Do not put the probe inside the cone or allow the cone to have direct contact with the probe.
- 2. Select the kVp, mA(s), and time normally used for an average adult patient thickness in routine bitewing radiography at that facility.
- 3. Make an exposure and document the radiation output in millirem. This direct measurement is the in-air exposure.



EXPOSURE REPRODUCIBILITY

$$C = \frac{s}{\overline{x}} = \frac{1}{\overline{x}} \left[\frac{\sum_{i=1}^{n} (x_i - \overline{x})^2}{n-1} \right]^{1/2}$$

Where:

s = Estimated standard deviation of the population.

 \overline{x} = Mean value of observations in sample.

 $X_i = i^{th}$ observation in sample.

n = Number of observations in sample.

Example:

The four (n) exposures (X_i) measured 409 mR, 387 mR, 391 mR, and 410 mR

STEP 1 Determine the mean value (\bar{x}) of the four exposures taken

(409 mR + 387 mR + 391 mR + 410 mR) ÷ 4 = 399.25 mR

STEP 2 Find the difference between each exposure and the mean value (\bar{x}) (disregard sign)

409.00 mR	387.00 mR	391.00 mR	410.00 mR
<u>-399.25 mR</u>	<u>-399.25 mR</u>	<u>-399.25 mR</u>	<u>-399.25 mR</u>
9.75 mR	12.25 mR	8.25 mR	10.75 mR

STEP 3 Square each of the differences

(9.75) ² =95.06	$(12.25)^2 = 150.06$
(10.75) ²⁼ 115.56	(8.25) ² = 68.06

STEP 4 Divide each number by 3 (n-1) and add the results

 $95.06 \div 3 = 31.69$ $150.06 \div 3 = 50.02$ $68.06 \div 3 = 22.69$ $115.56 \div 3 = \underline{38.52}$ 142.92

STEP 5 For s, determine the square root of the above number $\sqrt{142.92} = 11.95$

STEP 6 Divide s by the mean value (\overline{x})

11.95 \div 399.25 = .0299 = the coefficient of variation (C)

STEP 7 If $C \le 0.05$, the exposures are considered to be reproducible