

NEBRASKA DEPARTMENT OF HEALTH AND HUMAN SERVICES			
DIVISION OF PUBLIC HEALTH			
X-RAY PROGRAM			

### DENTAL EQUIPMENT PERFORMANCE EVALUATION

Manufacturer:		Serial Nur	nber:	
Coefficient of Variation:				
Measured:	mR	mR	mR	mR
Technical Factors Sele	ected: kVp _	mA / mAs	seconds	/ milliseconds / pulses
EXPOSURE REPRODUCIBILITYPASSFAILWhen all technique factors are held constant, the coefficient of variation of exposures for both manual and automatic exposure control systems must not exceed 0.05. This requirement applies to clinically used techniques.S = estimated standard deviation of the population X = mean value of observations in sample Xi= ith observation in sample n = number of observations in sample.				
Deviation:	%	%	%	%
Measured:				
Time Used for Testing:	·		s 🗌 milliseconds 🗌 pu	ulses
Manufacturer's Specifie	cations:	0	R ±10%	
TIMER The accuracy of the tin not obtainable, the time	ner must meet the mar er accuracy must be ±'	nufacturer's specificatic 10% of the indicated tir	<b>PASS</b> ons. If the manufacture ne with the testing per	<b>FAIL</b> r specifications are formed at 0.5 second.
Serial Number:		Location:		
Manufacturer:		Model Nur	nber:	
Control Panel Informat	ion			
Ion Chamber:	hin A Housing	External Probe		
Survey Instrument:	Ca	libration Date:		
Service Company:	Re	gistration Number: R	S	
Registrant:	Re	gistration Number:		
Evaluation Date:				

NEBRASKA			
Good Life. Great Mission.			
DEPT. OF HEALTH AND HUMAN SERVICES			

<b>KILOVOLT PEAK PASS FAIL</b> If the registrant possesses documentation of the appropriate manufacturer's kilovolt peak specifications, the radiation machine must meet those specifications. If the registrant does not possess documentation of the appropriate manufacturer's kilovolt peak specifications, the indicated kilovolt peak must be accurate to within ±10% of the indicated setting(s). For dental radiation generating equipment with fewer than three fixed kilovolt peak settings, the radiation machine will be checked at those settings.				
Manufacturer's Specification	ons:		DR ±10%	
Kilovolt Peak Used for Tes	ting:			
Measured:				
Deviation:	%	%	%	%
TUBE STABILITY PASS FAIL   The tube must remain physically stable during exposures. In cases where tubes are designed to move during exposure, the registrant must assure proper and free movement of the unit. FAIL				
COLLIMATION Field limitation must meet	the requirements o	of 180 NAC 21-007.07	<b>PASS</b> or 21-007.08	🗌 FAIL
Intraoral: Source-to-Skin Distance ( <i>AND</i> Beam limited to a diamete	SSD): cm r of 7 cm at minimu	Greater tl um SSD:	nan 18 cm: 🏾 🗌 Yes ] No	🗌 No
Panoramic:   Transverse axis: x-ray beam restricted to 0.0 inches of the imaging slit:   Yes   No   AND   Vertical axis: x-ray beam restricted to no more than 0.5 inches larger than the imaging slit:   Yes   No				
Cephalometric:   Rectangular collimation: x-ray field does not exceed 2.0% of the source-to-image receptor distance for the length or the width of the image receptor:   Yes   No   OR   Circular or polygon collimation: x-ray field does not exceed 2.0% of the source-to-image distance for the diagonal of the image receptor:   Yes   No				
For use on dental intraoral	systems on an av	erage adult patient thi	ckness in routine bitewi	ing radiography
Technique Factors:	kVp	mA / mAs	seconds / mill	iseconds / pulses
For Calculated Measurements Source to Skin Distance (S	ent only SSD):	cm Source to	Detector Distance (SD	)D):cm
Measured:r	nR 🗌 C	Calculated Measureme	ent Direc	t Measurement
Surveyor Name:		Surveyor Sig	gnature:	

# DETERMINING IN-AIR EXPOSURE MEASUREMENT FOR INTRAORAL DENTAL EXAMINATIONS

# A. CALCULATION

Note: Ion chambers may be located within the instrument housing rather than within an external probe. In this situation, the distance from the top surface of the housing to the ion chamber below must be known. If this type of instrument is used for the measurements, the inverse square law must be utilized for accurate results.

## IAE = Measured X (SDD ÷ SSD)<sup>2</sup>

Where:

IAE = in-air exposure Measured = indicated exposure on measuring instrument SDD = source (target) to detector (ion chamber) distance in centimeters SSD = source (target) to skin distance in centimeters

- 1. Place the tip of the cone within  $\frac{1}{2}$  inch from the housing of the measuring instrument.
- 2. Measure the distance from the source to the entrance/tube side surface of the housing. This is the SSD.
- 3. Determine the distance from the source to the ion chamber within the housing. This is the SDD.
- 4. Select the kVp, mA(s), and time normally used for an average adult patient thickness in routine bitewing radiography at that facility.
- 5. Make an exposure and document the radiation output in millirem.
- 6. Using the above formula, calculate the in-air exposure.

#### **B. DIRECT MEASUREMENT**

Note: Use this procedure only if an external probe (ion chamber) is available for the measurements.

- 1. Position the tube so the end of the cone is not greater than ½ inch from the probe. Do not put the probe inside the cone or allow the cone to have direct contact with the probe.
- 2. Select the kVp, mA(s), and time normally used for an average adult patient thickness in routine bitewing radiography at that facility.
- 3. Make an exposure and document the radiation output in millirem. This direct measurement is the in-air exposure.

# **EXPOSURE REPRODUCIBILITY**

$$C = \frac{s}{\overline{x}} = \frac{1}{\overline{x}} \left[ \frac{\sum_{i=1}^{n} (x_i - \overline{x})^2}{n-1} \right]^{1/2}$$

Where:

s = Estimated standard deviation of the population.

 $\overline{x}$  = Mean value of observations in sample.

 $X_i = i^{th}$  observation in sample.

n = Number of observations in sample.

Example:

The four (n) exposures (X<sub>i</sub>) measured 409 mR, 387 mR, 391 mR, and 410 mR

STEP 1 Determine the mean value (  $\overline{x}$  ) of the four exposures taken

(409 mR + 387 mR + 391 mR + 410 mR) ÷ 4 = 399.25 mR

STEP 2 Find the difference between each exposure and the mean value ( $\bar{x}$ ) (disregard sign)

409.00 mR	387.00 mR	391.00 mR	410.00 mR
<u>-399.25 mR</u>	<u>-399.25 mR</u>	<u>-399.25 mR</u>	<u>-399.25 mR</u>
9.75 mR	12.25 mR	8.25 mR	10.75 mR

STEP 3 Square each of the differences

(9.75) <sup>2</sup> =95.06	$(12.25)^2 = 150.06$
(10.75) <sup>2=</sup> 115.56	$(8.25)^2 = 68.06$

STEP 4 Divide each number by 3 (n-1) and add the results

 $95.06 \div 3 = 31.69$   $150.06 \div 3 = 50.02$   $68.06 \div 3 = 22.69$   $115.56 \div 3 = \underline{38.52}$ 142.92

STEP 5 For s, determine the square root of the above number  $\sqrt{142.92} = 11.95$ 

STEP 6 Divide s by the mean value (  $\overline{x}$  )

11.95  $\div$  399.25 = .0299 = the coefficient of variation (C)

STEP 7 If  $C \le 0.05$ , the exposures are considered to be reproducible