

Peter Iwen, PhD, D (ABMM), Director, Public Health Laboratory
Thomas J. Safranek, M.D., Nebraska State Epidemiologist
Samuel J. Pirruccello, M.D., Director, Client Services

981180 Nebraska Medical Center
Omaha, NE 68198-1180
Phone: (402) 559-2440
Toll Free: 866-290-1406
Fax: (402) 559-9497

Influenza Surveillance Specimen Collection and Transport Requirements

Web Link: <http://dhhs.ne.gov/publichealth/Documents/SpecimenCollectionRequirements.pdf>

Specimen: Nasopharyngeal washings or nasopharyngeal swab
Standard, contact, and airborne precautions are recommended for patient management; this includes collection of respiratory specimens. Nasopharyngeal washings are collected with normal saline, followed by a swab of the nasopharyngeal mucosa and both samples placed in the same viral transport media container. (See collection procedure below)

Volume: 2-3 mL nasopharyngeal washings and nasopharyngeal swab place in the same viral transport media container. Viral transport media container should be labeled with the patient first and last name, collection date and time, and specimen source.


Transportation: Follow enclosed FedEx packaging & shipping instruction. The specimens must be refrigerated from time of collection to arrival at NPHL, a pre-frozen gel pack is recommended. Samples must be received in lab within 72 hours of collection.

Courier service is provided by your local reference laboratory as identified on the test requisition. Your reference laboratory has agreed to participate in the Nebraska Public Health Partner Courier Network. Please contact NPHL Client Services at 402.559.2440 or 866.290.1406 if you are not sure of the specimen transport option available to your facility.

Specimen Collection Procedures

Nasopharyngeal Swab:

1. Using the Universal Viral Transport kit provided, insert the flocked flexible minitip through nares, parallel to the palate (NOT upwards) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient indicating contact with the nasopharynx. See Figure 1.
2. Gently, rub and roll the swab. Leave the swab in place for several seconds to absorb secretions before removing.
3. Place the swab back into viral transport media. NOTE: Combining the aspirate and the swab enhances viral recovery.
4. Transport refrigerated. Specimen must be received at NPHL within 72 hours of collection.
5. A video demonstrating the collection of a nasopharyngeal swab is available at: <http://www.youtube.com/watch?v=hXohAo1d6tk>.



Peter Iwen, PhD, D (ABMM), Director, Public Health Laboratory
Thomas J. Safranek, M.D., Nebraska State Epidemiologist
Samuel J. Pirruccello, M.D., Director, Client Services

981180 Nebraska Medical Center
Omaha, NE 68198-1180
Phone: (402) 559-2440
Toll Free: 866-290-1406
Fax: (402) 559-9497

Washing:

1. Cut off the distal end of the butterfly catheter (needle and butterfly) catheter extension set so that about 2-3 inches of tubing are left attached to the hub.
2. Draw up 2-3 mLs of sterile saline into a syringe.
3. Attach syringe to hub of butterfly catheter. Purge tubing with saline.
4. Put on gloves, gown, mask and eye protection.
5. Gently remove excess mucous from patient's nose. (If patient is an adult, ask the patient to gently blow nose. For pediatrics, a bulb syringe may be used to remove excess mucous.)
6. Position patient in supine position with the head of bed up 30°. The head should be turned to one side and tilted slightly backward.
7. Stabilize the patient's head and gently place the catheter into the nare. Placement should be in the nare (nasal wall). Depending on the size of the patient, this should be about 1-2 cm in adults and 0.5 cm to 1.0 cm in children (0.5 cm in neonates). See Figure 2.

8. Instill .5 - 2 mLs saline (.5 - 1 mLs for infants and children, 1 - 2 mLs for adults) into the nare and aspirate back mucous, saline and epithelial cells.
9. Repeat this process using the same syringe until the sample is cloudy or appears to hold cellular debris. (If the sample is inadequate, the process may be repeated on the opposite nare, using a second sterile syringe and tubing. Usually one nare is sufficient.)

NOTE: There may be some blood streaks in the mucous. This is normal and patients/parents should be told this is expected and will stop in a few minutes.

10. Transfer contents of tubing and syringe into transport media using the following process: Depress syringe plunger and express fluid from syringe and tubing into transport media. Then withdraw media/fluid back into syringe and tubing. Depress syringe plunger again, expressing fluid from syringe and tubing back into transport media.

NOTE: This is necessary to recover any cells or virus adhering to the tubing or syringe.

11. Transport refrigerated. Specimen must be received at NPHL within 72 hours of collection.




Figure 2



Peter Iwen, PhD, D (ABMM), Director, Public Health Laboratory
Thomas J. Safranek, M.D., Nebraska State Epidemiologist
Samuel J. Pirruccello, M.D., Director, Client Services

NEBRASKA

Good Life. Great Mission.

DEPT. OF HEALTH AND HUMAN SERVICES

981180 Nebraska Medical Center
Omaha, NE 68198-1180
Phone: (402) 559-2440
Toll Free: 866-290-1406
Fax: (402) 559-9497

2017-18 NPHL INFLUENZA SURVEILLANCE PROJECT SPECIMEN SHIPPING INSTRUCTIONS

Using Nebraska Public Health Partner Courier Network

1. Place viral transport media in biohazard bag and secure closure.
2. Make a copy of the completed NPHL Special Influenza Microbiology Requisition and keep for your records.
3. Fold the original completed NPHL Special Influenza Microbiology Requisition and place in outside plastic sleeve of the biohazard bag.
4. Place the “NPHL” and “ship refrigerated” sticker on the outside of the plastic sleeve.
5. Place entire package in your laboratory’s refrigerated storage area for usual courier pickup.
6. Courier to NPHL through your reference laboratory courier. The specimen may also be shipped with the NPHL or Regional Pathology Services courier if this service is provided to your facility.
7. Call NPHL client services with questions or concerns at 402.559.2440 or toll free at 866.290.1406.



Peter Iwen, PhD, D (ABMM), Director, Public Health Laboratory
Thomas J. Safranek, M.D., Nebraska State Epidemiologist
Samuel J. Pirruccello, M.D., Director, Client Services

NEBRASKA

Good Life. Great Mission.

DEPT. OF HEALTH AND HUMAN SERVICES

981180 Nebraska Medical Center
Omaha, NE 68198-1180
Phone: (402) 559-2440
Toll Free: 866-290-1406
Fax: (402) 559-9497

2017-18 NPHL INFLUENZA SURVEILLANCE PROJECT SPECIMEN SHIPPING INSTRUCTIONS

Using FedEx

Note: Personnel shipping specimens via FedEx should be familiar with packing and shipping requirements for Category B, Biological Substances.

1. Use only the Saf-T-Pak 302 “shipping system” including biohazard bag, absorbent, Tyvek Envelope and pre-labeled insulated shipper. These MUST be used together to be UN drop ship certified. Please keep the attached Data Sheet, packaging instructions, FedEx airbill and copy of NPHL Special Influenza Requisition with your procedures for 2 years, as this information may be requested by a CAP, FAA or DOT inspector.
2. Tighten lid on transport container to prevent leakage. Place viral transport media in biohazard bag with absorbent and secure closure.
3. Fold the original completed NPHL Special Influenza Microbiology Requisition and place in outside plastic sleeve of the biohazard bag.
4. Place biohazard bag into white Tyvek envelope and secure closure. Write list of contents on outside of envelope.
5. Place pre-frozen gel pack and white Tyvek envelope into insulated shipper.
6. Fit styrofoam lid on top of chest.
7. Seal shipping container with packing tape.
8. Use partially completed FedEx airbill (with NPHL shipping address & account number), fill out your facility’s shipping address and contact phone number.
9. Give to FedEx courier that same day.
10. Call NPHL client services with questions or concerns at 402.559.2440 or toll free at 866.290.1406.